







# Packaging and shipment of samples

TYPE OF SAMPLE	BLOOD SERUM	FAECES	SWAB	SURFACE SAMPLE	ORGANS	
<b>TEMPERATURE</b>	<b>ROOM TEMPERATURE</b>				<b>ROOM TEMPERATURE OR AT 4 °C</b>	<b>-20 °C TO -80 °C</b>
Primary receptacle: leakage proof container						
Packaging code P650 of ADR <sup>1</sup>						
Secondary package						
<b>UN3373</b>						
<b>BIOLOGICAL SUBSTANCE CATEGORY B</b>						
Outer package for transport						

<sup>1</sup> Accord Européen relatif au transport international des marchandises Dangereuses par Route.

Shipment of diagnostic samples can be offered through our DHL account. Please contact our Customer Service Center for assistance. Contact us by phone: +31 (0)24 361 5433 or by e-mail: [qmdiagnosics@qmdiagnosics.org](mailto:qmdiagnosics@qmdiagnosics.org)