















# Packaging and shipment of samples

TYPE OF SAMPLE		BLOOD SERUM	FAECES	SWAB	SURFACE SAMPLE	ORGANS	
TEMPERATURE		ROOM TEMPERATURE				ROOM TEMPERATURE OR AT 4 °C	-20 °C TO -80 °C
Primary receptacle: leakage proof container	Packaging code P650 of ADR <sup>1</sup>	 		 		  	
Secondary package	<b>ADR</b>  <b>UN3373</b>  <b>BIOLOGICAL SUBSTANCE CATEGORY B</b>						
Outer package for transport		 					

<sup>1</sup> Accord Européen relatif au transport international des marchandises Dangereuses par Route. Shipment of diagnostic samples can be offered through a DHL account. Please contact our Customer Service Center for assistance. Contact us by phone: +31 (0) 24 3615433 or by e-mail: [qmdiagnosics@qmdiagnosics.org](mailto:qmdiagnosics@qmdiagnosics.org)