

Packaging and shipment of samples

TYPE OF SAMPLE		BLOOD SERUM	FAECES	SWAB	SURFACE SAMPLE	ORGANS	
TEMPERATURE		ROOM TEMPERATURE				ROOM TEMPERATURE OR AT 4 °C	-20 °C TO -80 °C
Primary receptacle: leakage proof container	Packaging code P650 of ADR ¹	 		 		  	
Secondary package	ADR UN3373 BIOLOGICAL SUBSTANCE CATEGORY B						
Outer package for transport			 				

¹ Accord Européen relatif au transport international des marchandises Dangereuses par Route.

We offer shipment of diagnostic samples. Please contact our Customer Service Center for assistance. Contact us by phone: +31 (0)24 361 5433 or by e-mail: qmdiagnosics@qmdiagnosics.org